


# F&BFOCUS

SERVICING THE FOOD & BEVERAGE MARKET

Volume 1 • Spring 2015

**Advancing  
Microbiology  
Applications for  
Food & Beverage**

**Page 8** Solutions for Automated Water Pathogen Detection | **Page 16** Detection of Pathogenic *Vibrio parahaemolyticus* in Marine Food Using Lateral Flow Technology | **Page 20** Analysis of Free Chlorine in Poultry Chiller Water using DPD-Free Chlorine Reagent



## Microbiology Applications for Food & Beverage

Safety of food has always been a priority for food and beverage manufacturers. With the inception of FSMA the need for rapid, reliable testing has only increased. VWR is here to help with a broad array of media, rapid tests, consumables, and instruments to support all your needs.

Featured Inside	Page
Vacuum Filtration Systems.....	3
Environmental Sampling Swabs.....	3
Blender Sample Bags.....	4
Valor Food Scales.....	4
Paddle Blender.....	5
Class II Biosafety Cabinets.....	5
Nitrile 200 Powder-Free Gloves.....	6
Light-Duty Tissue Wipers.....	6
Microbiological Incubators.....	6
Food Sample Preparation Essentials.....	7
<b>Article</b> - Solutions for Automated Water Pathogen Detection.....	8-9
QC Microorganisms.....	10
Sterile Laboratory Water.....	10
Minimal Premixed E. coli Media.....	10
Pipettes and Bottletop Dispensers.....	11
<b>Pull-Out Supplement</b>	
pH Indicator Strips - Two Free Packs.....	12
pH Electrodes - Special Pricing.....	12
Microbiology Media - Free Samples.....	12
Metal Detectable Sampling Tools - Free Samples.....	12
<b>Pull-Out Supplement</b>	
Pipettor Carousel Stand - Free with Pipettor Purchase.....	13
Weighing Accessories - Free with Food Scale Purchase.....	13
VWR Custom Manufacturing Services.....	13
Ready-to-Use Media for Bacterial Counts.....	14
pH Electrodes for Food & Beverage Applications.....	15
<b>Article</b> - Detection of Pathogenic <i>Vibrio parahaemolyticus</i> in Marine Food Using Lateral Flow Technology.....	16-18
Lateral Flow Tests for Food Pathogen Detection.....	19
<b>Article</b> - Analysis of Free Chlorine in Poultry Chiller Water using DPD-Free Chlorine Reagent.....	20
Portable Spectrophotometer.....	21
HQd Water Quality Meters & Probes.....	21
Automatic Titrators.....	21
Microbiology Pump, Membrane Dispenser, & Filter Funnels.....	22
<b>VWRCATALYST Services</b> .....	22
Anaerobic Chambers.....	23
Metal Detectable Sampling Tools.....	23
Automated Water Pathogen Detection.....	24

**VWR has All You Need** for the food and beverage industry, offering the most extensive product portfolio and services supporting food analytical labs and production facilities.

**"All You Need" for Food Analysis Catalog:** comprehensive offering of products to support food testing

**VWR Custom Manufacturing Services:** providing high-quality chemicals and services for food quality testing

**VWR CATALYST™ Services Portfolio:** deliver process improvements to your lab or production facility

- **VWR Furniture catalog:** from design to project completion
- **Industrial Supplies catalog:** janitorial, MRO, and much more
- **Production Insight:** products and articles to support production
- **VWR Safety Diamond:** products to meet your safety needs



**VWR**  
COLLECTION



## VWR® Vacuum Filtration Systems featuring Pall Membranes

Designed for the sterile preparation of buffers, tissue culture media, microbiological media, and other biological fluids. Systems feature large diameter Pall membranes for guaranteed performance and maximum throughput.

Includes a filter funnel, 15 mL or 50 mL centrifuge tube with polypropylene screw cap, and vacuum hose adapter. Designed to filter fluids directly into the tube, eliminating the potential for contamination and spills during post-sterilization transfer.



Size	Pore Size, µm	Cat. No.	Case of 12
250 mL	0.1	89220-694	125.10
500 mL	0.1	89220-696	159.06
1 L	0.1	89220-698	200.20
250 mL	0.2	97066-200	100.37
500 mL	0.2	97066-202	148.09
1 L	0.2	97066-204	190.35
250 mL	0.45	97066-206	100.37
500 mL	0.45	97066-208	148.09
1 L	0.45	97066-210	190.35

Visit [vwr.com](http://vwr.com) to find full assembly with centrifuge tube options, as well as those with funnel only, bottle-top or centrifuge tube.

## EnviroMax® and EnviroMax Plus®

### The Ultimate in Environmental Sampling from the Swab Experts at Puritan

We spent time in the marketplace with different end-users, from food production to cosmetics to pharmaceutical firms. We saw how you work - listened to what you need - then went to our R&D team to design best-in-class environmental swabs.

- Oversized - perfect for sampling surfaces or large equipment parts
- Quick-turn open and leak-resistant cap
- Swabs are securely attached to caps
- Rigid paddle inside foam is ideal for reaching into crevices
- Maximized collection and elution
- EnviroMax Plus Swabs are pre-moistened with 50% neutralizing buffer & 50% 0.1% peptone water

#### Specifications

<b>Handle Length</b>	147.32 mm
<b>Overall Length</b>	152.09 mm
<b>Tip Diameter</b>	15.54 mm
<b>Tip Length</b>	56.38 mm



Description, Tip	Sterility	Swab Length, mm	Cat. No.	Price
<b>EnviroMax Dry Swabs</b>				
Pointed	No	147.32	89221-726	Pk. 25/ 54.04
Pointed	No	147.32	89221-728	Cs. 250/ 492.80
Pointed	Yes	147.32	89221-730	Pk. 25/ 58.49
Pointed	Yes	147.32	89221-732	Cs. 250/ 516.72
Rounded	No	145.03	89221-734	Pk. 25/ 54.04
Rounded	No	145.03	89221-736	Cs. 250/ 492.80
Rounded	Yes	145.03	89221-738	Pk. 25/ 58.49
Rounded	Yes	145.03	89221-740	Cs. 250/ 516.72
<b>EnviroMax Plus Pre-Moistened Swabs</b>				
Pointed	Yes	147.32	89221-742	Pk. 25/ 58.78
Pointed	Yes	147.32	89221-744	Cs. 250/ 508.39
Rounded	Yes	145.03	89221-746	Pk. 25/ 58.78
Rounded	Yes	145.03	89221-748	Cs. 250/ 508.39

# SAMPLE PREPARATION



## VWR® Blender Sample Bags

Standard bags are made of highly transparent polypropylene and have wide-sealed bottoms to prevent leaks and a tear-off sterile barrier top to ensure sterility up to the time of use.

FILTRA-BAG® bags are made of a mix of polyethylene and nylon. The dividing filter membrane allows for optimal solution and bacterial flow between the compartments during blending. Residue and semi-solid or solid substances remain in the compartment where the sample was inserted. Liquid can be extracted from the other compartment for an accurate aliquot that is free of solid particles.

SECURE-T® bags are made of heat extruded polyethylene tubing for internal sterility and to eliminate the need for side seals. Bags are malleable yet durable, with superior wall strength. May be used in general purpose blending or analysis testing.



Description	Capacity	Thickness, mil.	Cat. No.	Pack of	Case of
Transparent Polypropylene	650 mL	3	82007-698	500/ 133.78	1000/ 231.84
FILTRA-BAG, Open Top	1.63 L	3	10048-880	100/ 110.13	400/ 440.52
FILTRA-BAG, Safety Tab, Flat Wire Closure	710 mL	3	89085-570	100/ 82.25	400/ 284.99
FILTRA-BAG, Tear-Off Top	1.63 L	3	89085-572	100/ 91.95	400/ 317.96
SECURE-T, Tear-Off Top	1.63 L	4	11216-902	250/ 86.41	1000/ 299.42
SECURE-T, Tear-Off Top, Write-On Strip	1.63 L	3	89003-820	250/ 67.75	1000/ 258.95



## Meet the Newest Valor Food Scales

The new OHAUS Valor Series of compact food scales are reliable and durable products that will increase your food weighing and processing productivity. Dual displays, user-friendly menu and setup, along with the speed and accuracy of the results provide for a seamless and hassle-free user experience that is second to none.

**Buy a Valor 7000 and receive a Free Accessory!**  
(Choose Ethernet Kit, RS232 Kit, USB Kit, or In-Use Cover Set - See page 13 for more information and redemption details.)



Valor 2000

Valor 7000

OHAUS Model	Capacity, kg (lb)	Readability, g (lb)	Cat. No.	Each
<b>Valor 2000 Models</b>				
V22PWE1501T	1.5 (3)	0.0002 (0.0005)	10123-886	440.04
V22PWE3T	3 (6)	0.0005 (0.001)	10123-888	440.04
V22PWE6T	6 (15)	0.001 (0.002)	10123-890	440.04
V22PWE15T	15 (30)	0.002 (0.005)	10123-878	440.04
<b>Valor 4000 Models</b>				
V41XWE1501T*	1.5 (3)	0.0002 (0.0005)	10123-926	482.10
V41XWE3T*	3 (6)	0.0005 (0.001)	10123-896	482.10
V41XWE6T*	6 (15)	0.001 (0.002)	10123-898	482.10
V41XWE15T*	15 (30)	0.002 (0.005)	10123-900	482.10
<b>Valor 7000 Models</b>				
V71P1502T	1.5 (3)	0.05 (0.0001)	89497-046	551.14
V71P3T	3 (6)	0.1 (0.0002)	89497-048	551.14
V71P6T	6 (15)	0.2 (0.0005)	89497-050	551.14
V71P15T	15 (30)	0.5 (0.001)	89497-052	551.14
V71P30T	30 (60)	1 (0.002)	89497-054	551.14

\*Includes Stainless Steel Housing





## Paddle Blender

The Boekel Scientific Triple Mix Paddle Blender is designed for use in food safety testing and macro sample preparation applications. The unit is designed to be used with standard and filter-type food sample preparation bags. The unique three paddle design provides excellent homogenization and mixing.

The large process viewing window and easy-to-use digital control system allows users to view the sample status without opening the door. The unit features an easy to clean design with removable paddles and drip tray. To facilitate use in high throughput labs, an Auto-Run feature is included, allowing users to set and repeat the same times and rpm's in a continuous fashion.

- Three paddle design for homogenization and mixing
- Easy to use digital control system with an auto-run feature
- Large process viewing window
- Easy to clean design with drip tray and removable paddles
- Two-year Warranty

### Specifications

<b>Speed Range</b>	75 – 300 RPM; 5 RPM increments
<b>Capacity</b>	80 – 400 mL bags
<b>Timer</b>	0 – 999 minutes 59 seconds
<b>Operating Temperature Range</b>	10°C to 30°C

Electrical	Cat. No.	Each
115V, 50/60Hz, 150W	10040-414	5529.94
230V, 50/60Hz, 150W	10040-416	5781.93



## Rediship Purifier® Logic® + Class II Biosafety Cabinets Designed for the HUMAN Experience

These Biosafety Cabinets raise the bar for not only safety, but comfort, too. They were designed for the human body and our habits. We call it Inclination™ Technology — features that reinforce what humans are inclined to do. We believe the more comfortable you are, the safer you will be. Features include:

**Performance:** Constant Airflow Profile™ (CAP) technology works with the ECM for precise airflow control

**Safety:** MyLogic™ OS visible on a large, color LCD display; NSF-listed

**Comfort:** Features include true line-of-sight display, curved inlet grille for forearm support, and flush-mounted electrical receptacle covers with dampened hinges for easy opening and closing

**Sustainable Design:** Contributes toward a building's energy use reduction LEED points; has 94.9% recyclable content when the cabinet's lifetime is completed

Description	Dimensions, W x D x H, in.	Cat. No.	Each
3' with Stand, 10" Sash	42-5/16 x 36-1/2 x 13-5/8	89413-124	10,185.00
4' with Stand, 10" Sash	54-5/16 x 48-1/2 x 19-7/8	89413-128	11,475.00
5' with Stand, 10" Sash	66-5/16 x 60-1/2 x 26-7/8	89413-132	12,780.00
6' with Stand, 10" Sash	78-5/16 x 72-1/2 x 32-7/8	89413-136	14,140.00

Contact your VWR Sales Representative for more information on additional options and accessories. Purifier Logic+ and accessories are available for delivery through the VWR REDISHIP Program.



## VWR® Nitrile 200 Powder-Free Gloves

Manufactured from 100% premium nitrile, VWR Nitrile 200 Powder-Free Examination Gloves provide hand protection against biohazards and a broad range of chemicals. 100% latex-free, these gloves eliminate adverse reactions associated with natural rubber latex proteins. Gloves are rigorously tested and manufactured under stringent process controls according to ISO 9001 and ISO 13485.



### Specifications

Length	23 cm (9")
Cuff Thickness	2.4 mil
Finger Thickness	3.9 mil
Palm Thickness	2.8 mil

Size	Cat. No.	Box of 200	Case of 10
X-Small	89428-746	31.09	242.41
Small	89428-748	31.09	242.41
Medium	89428-750	31.09	242.41
Large	89428-752	31.09	242.41
X-Large	89428-754	31.09	242.41
XX-Large	89428-756	31.09	242.41



## VWR® Light-Duty Tissue Wipers

- Virgin fiber tissue
- Absorbent, strong, and low-linting
- Antistatic, film-windowed boxes

Wipers come in antistatic film-windowed dispenser boxes. Custom-size dispensers hold wiper box securely, providing one-handed access to wipers. Dispensers are made from 6mm (1/4") white PVC and can be wall-mounted or placed on a bench. Double-faced tape allows semi-permanent installation and holds dispenser securely. Medium and large dispensers feature keyhole openings for screws, providing for permanent installation. Item **82003-820** features perforated top and bottom access for all standard dispensers.



Custom-sized dispensers hold wiper box securely for one-handed access

Packaging	W x L, cm	Thickness	Cat. No.	Pack of	Case of
280/Box	11.4 x 21.3	1-Ply	82003-820	280/ 5.79	16,800/ 270.30
140/Box	35.6 x 42.4	1-Ply	82003-822	140/ 18.86	2100/ 187.97
90/Box	35.6 x 42.4	3-Ply	82003-824	90/ 28.09	1350/ 288.03

Dispensers	Cat. No.	Each
Small, for 11.4 x 21.3L cm Wipers	14233-756	14.76
Medium, for 35.6 x 42.4 cm Wipers	14233-758	22.04
Large, for 35.6 x 42.4 cm Wipers	14233-760	27.19



## VWR® Microbiological Incubators

VWR Forced Air incubators provide an even higher temperature uniformity and stability for reproducible results in microbiology and research labs. Compared to gravity convection units the temperature recovers more quickly after opening door and when replacing samples. Optimal conditions are maintained even with larger loads.

VWR Gravity Convection incubators with unique airflow design provide best conditions for delicate samples. Drying out is kept at a minimum level, even in open containers. Besides incubation in microbiology, bacterial and yeast research, the VWR incubators are used in food/beverage/pharmaceutical testing, and for heated storage.



- Two convection choices - gravity or forced air
- Three table-top sizes: 60L, 100L, 180L and a 400L floor model
- Intuitive, easy to use controller
- Excellent temperature uniformity and stability
- Stainless steel interior with round corners for easy cleaning
- Flexible shelving

Volume, cu. ft.	Uniformity at 37°C	Stability at 37°C	Shelves Supplied*/Max	Cat. No.	Each
<b>Forced Air Incubators</b>					
2.3	± 0.2°C	± 0.1°C	2/13	89511-424	2301.33
3.6	± 0.3°C	± 0.1°C	2/16	89511-426	3039.11
6.3	± 0.4°C	± 0.1°C	2/19	89511-428	3414.42
13.4, Free-standing	± 0.3°C	± 0.2°C	2/39	89511-430	6248.73
<b>Gravity Convection Incubators</b>					
2.6	± 0.6°C	± 0.2°C	2/13	89511-418	1798.51
4	± 0.6°C	± 0.2°C	2/16	89511-420	2317.28
6.85	± 0.6°C	± 0.2°C	2/19	89511-422	2617.05

\*Additional perforated shelves and wire shelves are available, visit [vwr.com](http://vwr.com).

# All You Need For Food Sample Preparation


Our range of consumables and equipment for sample preparation includes not only a comprehensive selection of products, but also a choice of materials to meet every need.

- ▶ Beakers
- ▶ Benchtop equipment
- ▶ Blender **NEW!**
- ▶ Blender bags
- ▶ Bottles
- ▶ Bag clips
- ▶ Conical and straight containers
- ▶ Centrifuge tubes
- ▶ Diluter **NEW!**
- ▶ Disposable glass tubes
- ▶ Erlenmeyer flasks
- ▶ Filters
- ▶ Loops
- ▶ Needles
- ▶ Petri dishes
- ▶ Pipets and pipet controllers



**GOSSELIN™**

*A Corning Brand*



**FALCON®**

*A Corning Brand*

**CORNING**

**FALCON**

**AXYGEN**

**GOSSELIN**

**PYREX**

**Warranty/Disclaimer:** Unless otherwise specified, all products are for research use only. Not intended for use in diagnostic or therapeutic procedures. Corning Life Sciences makes no claims regarding the performance of these products for clinical or diagnostic applications.

For a listing of trademarks, visit us at [www.corning.com/lifesciences/trademarks](http://www.corning.com/lifesciences/trademarks). All other trademarks are the property of their respective owners.

©2015 Corning Incorporated. All rights reserved.

## Solutions for Automated Water Pathogen Detection

by Dr. R. Stephen Brown, Associate Professor, Department of Chemistry and the School of Environmental Studies at Queen's University, Ontario, CA



### Safe Water is a Critical Need

Water plays a central role in our lives, and access to safe water is essential for everyone – from those living in the world's largest cities to the smallest towns and as a key element of industrial processes producing everything from the most basic foods to today's most advanced microelectronics. Water destined for human and industrial consumption is required to be quality tested for a wide variety of contaminants, including potentially harmful microorganisms. While many water quality parameters such as pH, chlorine and turbidity can be measured in near real-time by on-line measurement instrumentation, microbiological testing presents a unique challenge. The requirement to detect a single *E. coli* cell in a 100mL water sample has been compared to the challenge of finding a single coffee bean in 40,000 Olympic-sized swimming pools.


The majority of microbiological water quality tests rely on traditional microbiological methods that were developed decades ago. Water samples are sent to a laboratory where they are filtered and subsequently placed into a growth culture media that is typically incubated for 24 hours. After the incubation period, each sample is inspected by a trained technician for signs of bacterial growth. Chemical additives in the growth media cause colonies of specific target bacteria to change color or to appear fluorescent under ultraviolet light. While these methods are the current gold standard, they require 24-48 hours for samples to be transported, analyzed and the laboratory results to become available. In the event of contamination, delays in receiving test results increase the risk of public exposure, requiring costly remedial action or product recalls.

### The Innovation: Introducing the TECTA Automated Microbiology Solution

The impact of highly-publicized events involving contamination of municipal water supplies has driven the need to develop new microbiological test methods that are faster and easier to use than traditional methods, while maintaining the sensitivity and reliability that have been trusted for decades.

In 2001, a consortium of university researchers and water quality experts proposed an innovative approach – a novel way of automating the test by using a polymer-based optical sensor to detect the same types of fluorescent indicators of bacteria used in the trusted current methods. Pathogen Detection Systems, Inc. was subsequently formed to provide municipalities and industries worldwide with the automated microbiology solutions that were developed as a result of this research.





**U.S. EPA Approved**

The innovative capabilities of the TECTA system are made possible by a patented polymer based optical sensor, the “Optical Partition Element”, that is built into every test cartridge. Each cartridge contains pre-measured amounts of growth media that support the enrichment of any target bacteria that are present in the sample. As target bacteria such as *E. coli* or Coliforms begin to multiply, they emit a specific enzyme that interacts with a proprietary chemical substrate in the cartridge, releasing fluorescent molecules from the substrate. These fluorescent indicators rapidly move from the water sample into the polymer optical sensor located within the cartridge, enabling automated detection by a low-cost ultraviolet optical detection system that is built into the TECTA instrument.

The novel combination of the polymer-based optical sensor and ultraviolet optical detection system provides several key

advantages. The fluorescent indicators are extracted and concentrated within the polymer of the optical sensor, facilitating both rapid detection at the earliest possible time and eliminating the risk of off-color samples or turbidity within the water sample obscuring the fluorescent indicator. When placed in the TECTA instrument the test is monitored continuously throughout the incubation process – providing alerts of contaminated samples as soon as possible and eliminating the requirement to wait until the end of a fixed incubation process (typically 18-24 hours). This unique early alerting capability allows the system to provide results within 2-18 hours, depending on the level of sample contamination. In addition to a “presence/absence” result, the system is capable of providing an estimate of the number of bacteria that were present in the original sample – an important indicator of the level of severity of an adverse microbiological test result.

*Dr. R. Stephen Brown's research is focused on instrumentation and methodology related to water quality monitoring and environmental toxicology. He has published over 80 peer-reviewed papers, five book chapters, and has six issued patents. He has given many invited seminars and conference presentations on his research and on technology transfer. He co-founded the company Pathogen Detection Systems, now a subsidiary of Veolia Environment in the ENDETEC group. He continues to work with ENDETEC as Chief Scientist to develop and commercialize new water monitoring technologies.*



## TECTA™ Automated Microbiology System

TECTA is the first automated microbiological testing system approved by U.S. EPA for regulatory compliance testing of municipal drinking water systems required under the Total Coliform Rule.

The system saves precious time and money by being able to test onsite without the need for a microbiologist. It requires much less manual labor and use of disposable products – the test cartridge is pre-filled with all required test reagents eliminating the need for any handling, dilution, or mixing of reagents for test samples.

- Shows positive results with single cell sensitivity and provides an estimate of the number of bacteria present
- Simplifies the FSMA process, easily handling the increased number of tests
- Approved by the U.S. EPA for drinking water regulatory compliance monitoring

Description	Cat. No.	Price
TECTA Automated Rapid Microbial Detection System	10218-898	Ea./ 19,750.00
Combined <i>E. Coli</i> and Total Coliform Test, 100 mL	10220-374	Box 48/ 472.00



## Epover™ QC Microorganisms for Multiple Application Enumerated Challenges



Epover QC microorganisms can be used to perform many essential quality control functions including microbial detection and enumeration, equipment calibration, method validation, bioburden determination, antibacterial effectiveness and lethality testing. Available in concentrations ranging from  $10^2$  to  $10^8$  CFU per pellet, Epover lyophilized microorganism preparations are extremely versatile. Epover pellets can be used individually, or several pellets can be combined for a mixed microorganism population challenge.

- Easily manipulated to deliver a wide variety of CFU levels
- Packaged in a resealable vial containing ten pellets of a single quantitative QC microorganism strain
- Saves time and reduces labor by delivering a specific number of colony forming units (CFU)
- Doesn't require frozen conditions - lower shipping costs
- Documentation of phenotypic, macroscopic, and microscopic test results
- Peel-off Certificate of Assay included for quality control documentation; meets USP guidelines
- Lower price than competitive products


**Over 45 strains available!**

Visit [vwr.com](http://vwr.com) or talk to your VWR Sales Representative to learn more.



**Same Pure Product.  
New Pure Brand.**

## Introducing a NEW VWR Life Science Brand!



Together, with VWR's expertise in delivering solutions and AMRESCO's pure product, we bring you our manufactured specialized solutions to better enable your science under one name - VWR.

Sterile Water	Size	Cat. No.	Price	
Low Endotoxin, Tissue Culture Grade	500 mL	VWRL0200-0500	Ea./ 9.95	Cs. 10/ 84.50
Low Endotoxin, Tissue Culture Grade	1 L	VWRL0200-1000	Ea./ 17.95	Cs. 6/ 101.50
Purified, Molecular Biology Grade	500 mL	VWRL0201-0500	Ea./ 7.95	Cs. 10/ 67.50
Purified, Molecular Biology Grade	1 L	VWRL0201-1000	Ea./ 13.95	Cs. 6/ 81.50
Nuclease-Free	100 mL	97062-790	Ea./ 19.33	—
Nuclease-Free	500 mL	97062-794	Ea./ 34.21	—

## Minimal Premixed *E. coli* Media

- Save time using premixed and ready-to-use media
- Homogenous, high-quality powders

Culture of bacteria is an important area in genetic and molecular biology research. Researchers need a variety of microbiological growth media that will support organisms in an efficient and productive manner. Selecting the media that will best serve the necessary requirements will allow an increase in both the bacterial yield and DNA expression.

Description	Size, g	Cat. No.	Each
H Medium Broth	500	97063-460	76.80
M63 Medium Broth	100	97063-466	13.67
M63 Medium Broth	500	97063-468	47.50
M9CA Medium Broth	100	97063-396	12.23
Tryptone Broth	100	97063-422	17.64
Tryptone Broth	500	97063-424	75.59




# Complete Pipetting

## From pipettes to calibration to bottletop dispensing

The Eppendorf complete liquid handling solution offers your lab a single solution for superior quality tips, calibration services and liquid handling products.

- > Eppendorf liquid handling products meet the highest needs for precision, ergonomics and robustness
- > Eppendorf offers a variety of tips in the purity grade your lab needs
- > Eppendorf also offers calibration services for all Eppendorf and non-Eppendorf brand pipettes



Description	Volume Range	Cat. No.	Each
Research® Plus Pipette, Adjustable	1000–10,000 µL	89125-310	394.21
Reference® 2 Single Channel Pipettor, Adjustable Volume	0.25–2.5 mL	10032-684	433.70
Xplorer® Plus Electronic Pipette	0.5–10 mL	89234-006	1033.82
Repeater® M4 Mechanical Dispenser Starter Kit*	1 µL–10 mL	10054-496	580.83
Top Buret™ H Digital Bottle Top Dispenser	0.01–999.9 mL	89133-974	1221.84
Varispenser® Bottletop Dispenser	1–5 mL	89133-950	512.46
Varispenser Bottletop Dispenser	10–50 mL	89133-956	727.26

\*Starter Kit includes Combitip® Rack, Combitip® assortment pack 1 µL to 1 mL dispensing volume (for use with Combitip® advanced 0.1–50 mL).

# MONEY SAVING OFFERS

Visit [us.vwr.com/promotions](http://us.vwr.com/promotions) and enter promo code in the "FIND" box for promo details and redemption information. See promotion terms and disclaimer at [us.vwr.com/promotions](http://us.vwr.com/promotions).

## MColorpHast™ pH Indicator Test Strips

Promo Code  
4307

### Free pH Strips

- Convenient method for accurate colorimetric pH indication
- Covalently bound indicator dyes prevent bleeding even in strong alkaline solutions
- Can be immersed in samples for extended periods, so that even weakly buffered solutions can be accurately tested without contaminating the sample



Buy a six-pack of the following pH indicator strips

Description	pH Measuring Range	Cat. No.	Each
Universal Indicator	0-14	EM1.09535.0007	123.02
Intermediate Indicator	0-6.0	EM1.09531.0007	123.00
Intermediate Indicator	5.0-10.0	EM1.09533.0007	123.00
Special Indicator	4.0-7.0	EM1.09542.0007	123.00
Special Indicator	6.5-10.0	EM1.09543.0007	123.00

No. of Test Strips: 600 (6 packs of 100)

Get two packs of the following pH indicator strips

Description	pH Measuring Range	Cat. No.	Each
Universal Indicator	0-14	EM1.09535.0001	Free
Intermediate Indicator	0-6.0	EM1.09531.0001	Free
Intermediate Indicator	5.0-10.0	EM1.09533.0001	Free
Intermediate Indicator	7.5-14.0	EM1.09532.0001	Free
Intermediate Indicator	2.0-9.0	EM1.09584.0001	Free
Special Indicator	0-2.5	EM1.09540.0001	Free
Special Indicator	2.5-4.5	EM1.09541.0001	Free
Special Indicator	4.0-7.0	EM1.09542.0001	Free
Special Indicator	6.5-10.0	EM1.09543.0001	Free
Special Indicator	11.0-13.0	EM1.09545.0001	Free

No. of Test Strips: 100

## Thermo Scientific Orion ROSS Ultra pH Electrodes

Promo Code  
4323

### Save 30%

Consistent Confidence in Your pH Measurements

For a limited time, instantly save 30% on selected Thermo Scientific™ Orion™ ROSS Ultra™ pH electrodes.

Thermo Scientific Orion ROSS Ultra pH electrodes stabilize quickly, provide highly accurate and precise measurements, and have an extended sensor life.

See page 15 in this issue for more information.

Thermo Scientific Orion ROSS Ultra Electrodes	Cat. No.	Each	
Glass, Refillable pH/ACT Triode™	MP89187-024	460.07	Save 30%
Epoxy, Refillable pH/ACT Triode	MP97011-870	423.78	
Epoxy, Refillable pH/ATC Triode, 3 m Cable	MP97011-872	361.31	
Epoxy, Low-Maintenance Gel pH/ATC Triode, 1.5 m Cable	MP97015-738	474.37	
Epoxy, Low-Maintenance Gel pH/ATC Triode, 3 m Cable	MP97015-740	491.29	Special Introductory Pricing
Epoxy, Low-Maintenance Gel pH/ATC Triode	10324-578	284.40	
Epoxy, Refillable pH/ATC Triode	10324-580	332.40	
Glass, Refillable pH/ATC Triode	10324-314	334.80	

pH Range: 0 to 14 for all featured electrodes, and all have waterproof BNC connectors

## Compact Dry™ Microbiology Media

Promo Code  
4312



Compact Dry is a ready-to-use system reducing the time needed to perform microbial testing on food, beverage, meat, cosmetics, and raw materials. See page 14 for more information on this simple and effective system for colony counting.



Get a FREE sample sleeve of your choice - no purchase necessary!

See page 14 in this issue for more information.

**Free Sample!**

Description	Used For	Cat. No.	Pack of 100	Sample (4 Plates)
Compact Dry CF	Detection of Coliforms	10145-952	152.60	Free
Compact Dry EC	Detection of Coliforms and E. coli	10145-954	185.15	Free
Compact Dry TC	Total Viable Bacterial Count	10145-968	129.15	Free
Compact Dry XSA	Detection of Staphylococcus aureus	10145-970	217.39	Free
Compact Dry YM	Differentiation of Yeasts and Mold	10145-972	202.33	Free

## Bel-Art Sterileware Senseable Scoop

Promo Code  
4206

**FREE SAMPLE PACK of Sterileware® Metal Detectable Sampling Tools includes:**

- (1) 89233-982 Blue, 60 mL Sterileware Senseable Scoop;
- (1) 89093-480 Black, 125 mL Sterileware Senseable Scoop;
- (1) 89233-986 Blue, 250 mL Sterileware Senseable Scoop



FDA grade – safe for food, drug, and cosmetic contact. See page 23 for product ordering information.

Offer expires 8.31.2015

**Scienceware**  
BEL-ART PRODUCTS

**Free Sample Pack!**

## Free Carousel Stand

Promo Code  
4332

value \$148!

with Purchase of an Eppendorf Research Plus or Reference 2 Six-Pack

**Research plus**—Remarkably light both in terms of weight and pipetting forces setting new standards for ergonomic operation. It also features a spring loaded tip cone for secure tip attachment and low ejection forces and is fully autoclavable.



**Reference 2**—In addition to all the features above, the Reference 2's one-button design minimizes the risk of aerosol contamination and reduces thumb strain even further.

6-Pack includes six single-channel, adjustable volume pipettes, one full box of Eppendorf epT.I.P.S. for each pipette (excludes 5 and 10 mL pipettes) and one carousel stand.

Description	Cat. No.	Each
<b>Reference 2, 6-Pack</b>		
0.5–10, 10–100, 30–300, 100–1000 µL, 0.5–5 mL, and 1–10 mL	89131-967	2070.00
0.1–2.5, 0.5–10, 2–20, 10–100, 20–200, and 100–1000 µL	89126-241	2070.00
<b>Research plus, 6-Pack</b>		
0.5–10, 10–100, 30–300, 100–1000 µL, 0.5–5 mL, and 1–10 mL	89429-358	1800.00
0.1–2.5, 0.5–10, 2–20, 10–100, 20–200, and 100–1000 µL	89126-236	1800.00

Offer expires 6.30.2015.

## Valor<sup>®</sup> 7000 Food Weighing Scales

Promo Code  
4331



Buy an OHAUS Valor<sup>®</sup> 7000 Scale and Receive a FREE Accessory. See page 4 for qualifying catalog numbers and ordering information.



The Valor 7000 compact food scale is built with features designed to improve food weighing and portioning productivity.

- Dual displays and touchless sensors
- User-friendly menu and setup
- Fast and accurate results

up to \$250 value

Choose One FREE Accessory	Cat. No.	Each
Ethernet Kit	89497-036	109.28
RS232 Kit	89497-038	109.28
USB Kit	89497-040	109.28
In-Use-Cover, Set (10)	—	250.00

Must redeem within 6 weeks of purchase date; limit one per lab. Offer expires 12.31.2015, or while supplies last.

## VWR Provides Custom Manufacturing Services for Food Quality Testing

The Right Ingredients for Faster Speed to Market



VWR provides high-quality chemicals and services, customized to fit the unique needs of food quality testing. Our global footprint and experience in molecular and diagnostic applications will ensure **rapid speed to market** so that your products can get to where they are needed most.

### OUR CORE CAPABILITIES SUPPORTING FOOD QUALITY TESTING NEEDS INCLUDE:

- Aseptic manufacturing
- Precision small-volume filling
- Cleanroom manufacturing and filling suites
- Handling and filling of temperature- or light-sensitive materials
- Powder handling, including particle size reduction, blending, and filling
- Handling and segregation of biological and animal origin components
- Homogeneous blending of complex powder blends
- Handling of antibiotic powders and solutions
- ISO 13485 compliant sites

Focus on your core business, and let VWR take care of sourcing, production, testing, and documentation of your chemicals. For more information, call **1.800.932.5000** or visit [vwr.com/manufacturing](http://vwr.com/manufacturing).

### Global Manufacturing Sites

Seven VWR manufacturing sites provide a global footprint to accommodate life science applications in biotechnology, biopharmaceutical, diagnostics, and research, encompassing research and development to full-scale production.



## Compact Dry™ for Bacterial Counts

Hardy Diagnostics, an FDA licensed and ISO certified biomedical manufacturer, has developed Compact Dry, a simple method for cultivating and counting common food-borne microorganisms.

Compact Dry is ready-to-use chromogenic medium for performing coliform counts, *E. coli*, yeasts and molds, as well as total aerobic bacterial counts. It contains dehydrated media and a water-soluble gelling agent in a non-woven cloth matrix. The medium is instantly hydrated when inoculated with a sample, and capillary action diffuses the sample evenly over the matrix to form a gel within seconds.

**Long shelf life (18 months) and can be stored at room temperature**

A diluted food or beverage sample (1mL) is added to the center of the device, which diffuses automatically and evenly throughout the plate. No spreader is required. The device is then incubated and colonies grow with specific colors for easy organism differentiation and counting.


**Compact Dry CF** - easy tool for the detection of coliforms

**Compact Dry EC** - detects both *E. coli* and coliforms


**Compact Dry TC** - medium for total viable bacterial count

**Compact Dry XSA** - determines presence of *Staphylococcus aureus*

**Compact Dry YM** - differentiates yeasts and mold by color



Compact Dry™ CF  
(Coliforms). Cat. No.  
10145-952




Compact Dry™ EC  
(*E. coli* & Coliforms).  
Cat. No. 10145-954



Compact Dry™ TC  
(Total Count). Cat. No.  
10145-968



Compact Dry™ XSA  
(*S. aureus*). Cat. No.  
10145-970



Compact Dry™ YM  
(Yeast & Mold). Cat.  
No. 10145-972



Description	Incubation Temperature	Incubation Time	Cat. No.	Pack of 100
Compact Dry CF	35°C	18-24 hr.	10145-952	152.60
Compact Dry EC	35°C	20-24 hr.	10145-954	185.15
Compact Dry TC	35°C	44-48 hr.	10145-968	129.15
Compact Dry XSA	35°C	24 hr.	10145-970	217.39
Compact Dry YM	25°C	3-7 days	10145-972	202.33

Packaging: Four plates per mylar sleeve, 25 mylar sleeves per pack.

# Consistent Confidence in Your pH Measurements

**Thermo**  
SCIENTIFIC



When measuring pH in the food and beverage industry, an incorrect reading can cost you time and money. It is critical that the pH measurements you perform day-to-day are fast, accurate and reproducible. Thermo Scientific™ Orion™ ROSS Ultra™ Triode™ pH/ATC electrodes stabilize quickly, provide highly accurate and precise measurements, and have an extended sensor life.

Why take chances with lower quality?

Thermo Scientific Orion ROSS Ultra Electrodes	Body Material	Junction	Connectors (pH/Temperature)	Cat. No.	Each	
Refillable pH/ATC Triode	Glass	Ceramic	BNC, 8-pin MiniDIN	MP89187-024	460.07	<b>Instantly Save 30%</b>
Refillable pH/ATC Triode	Epoxy	Glass Fiber	BNC, 8-pin MiniDIN	MP97011-870	423.78	
Refillable pH/ATC Triode, 3 m Cable	Epoxy	Glass Fiber	BNC, 8-pin MiniDIN	MP97011-872	361.31	
Low-Maintenance Gel pH/ATC Triode, 1.5 m Cable	Epoxy	Glass Capillary	BNC, 8-pin MiniDIN	MP97015-738	474.37	
Low-Maintenance Gel pH/ATC Triode, 3 m Cable	Epoxy	Glass Capillary	BNC, 8-pin MiniDIN	MP97015-740	491.29	<b>Special Introductory Pricing</b>
Low-Maintenance Gel pH/ATC Triode	Epoxy	Glass Capillary	BNC, RCA	10324-578	284.40	
Refillable pH/ATC Triode	Epoxy	Glass Fiber	BNC, RCA	10324-580	332.40	
Refillable pH/ATC Triode	Glass	Ceramic	BNC, RCA	10324-314	334.80	

All featured electrodes have a pH range of 0 to 14, have waterproof BNC connectors, and are compatible with TRIS buffers.



**For a limited time, instantly save 30% on selected electrodes (items with special 'MP' catalog numbers). Special pricing expires 7.31.2015.**

# Detection of Pathogenic *Vibrio parahaemolyticus* in Marine Food using Lateral Flow Technology

Lisa John, Jörg Slaghuis, & Heike Wulff; Merck Millipore, BioMonitoring R&D, Darmstadt, Germany



## Introduction

Foodborne illnesses impact consumer safety on a global scale. It is estimated that one in three people in industrialized countries may be affected by foodborne illness per year. In the United States, 48 million people get sick, 128,000 are hospitalized, and 3,000 die from foodborne illness annually. In the EU, 320,000 cases of foodborne illness are documented each year, although the real number is likely much higher.

Traditional microbiological methods for detection of pathogens in food can require up to five days to obtain a simple yes/no result. This time-consuming process slows the workflow, holding the food in quarantine and preventing its release. This can result in a considerable delay before products can be put into the market. Immunoassays based on the principle of lateral flow technology allow for convenient detection of pathogens within 24-48 hours, depending on parameter. These tests are available for a broad range of pathogens and follow a simple “pregnancy test” design to provide results in a quick, readable format and deliver definite results in as little as twenty minutes after sample enrichment. Lateral flow tests offer all the benefits of traditional testing methods with the addition of simplicity, speed, reliability and convenience. When used as part of a monitoring program, they allow streamlining of testing protocols, ensuring the safety of finished products and shortening holding times.


This article describes development and evaluation of a lateral flow test for pathogenic *Vibrio parahaemolyticus*, a major cause of foodborne illness throughout the

world, primarily associated with consumption of contaminated raw or undercooked seafood. (Please note that the lateral flow test for *V. parahaemolyticus* is not commercially available.)

Approximately 4500 cases of *V. parahaemolyticus* infection are reported each year in the United States. Numbers are expected to increase worldwide due to greater consumption of raw seafood and the globalization of seafood trade.

The thermostable direct hemolysin (TDH) toxin is known as the major virulence factor of *V. parahaemolyticus*. Standard detection methods of *Vibrio parahaemolyticus* vary by country, but in all cases are labor-intensive and require three to seven days for results. Because raw seafood quickly experiences deterioration, rapid detection methods are necessary for effective identification of possible contamination.

For this application, a Gold Labelled ImmunoSorbent Assay (GLISA), an immuno-chromatographic rapid test based on lateral flow technology (**Figure 1**), was used. The lateral flow assay (LFA) detects the toxin TDH using monoclonal gold-labelled antibodies. If the antigen is present, it reacts with the gold-labelled toxin-specific antibodies and migrates to the binding zone. The gold-labelled toxin-specific antibodies then link to a second specific antibody. Due to the gold-labelling, a distinct red line is formed. The rest of the sample continues to migrate to the control zone and links to a third antibody-specific antibody. The red line formed in the control zone demonstrates that the test is functioning correctly.



**Figure 1:** Lateral flow tests are designed to screen for the presence of pathogenic and toxigenic bacteria in food matrices. The kits are based on immunochromatographic principles and use antibody-linked colloidal gold particles to react specifically with its complementary antigenic determinant to provide a visual reaction read-out.

## Methods

Three studies were performed to evaluate the assay for TDH toxin:

- Limit of detection
- Inclusivity/exclusivity
- Evaluation with artificially contaminated food samples

To establish the limit of detection, four different TDH positive strains of *V. parahaemolyticus* pure cultures were diluted and tested with the LFA. Pure TDH was also tested.

Inclusivity and exclusivity of the LFA were evaluated by testing a total of 102 isolates and reference strains. Bacteria were cultured in Peptone water (acc. to ISO 6579) plus 2% NaCl pH 8.5 or in CASO Broth for 18 - 24 hours at 37 °C. A total of 160 µl of suspension was transferred onto the sample port of the test device. The result was read after 30 minutes.

For evaluation with artificially contaminated food samples, fish and seafood products (oysters, shrimp and sushi (e.g. nigiri and maki with salmon and tuna), n=90 total) were spiked with a TDH-positive *V. parahaemolyticus* strain and analyzed comparatively by the developed


test and the reference method according to ISO/TS 21872-1:2007 (**Figure 2**). For inoculation, stressed and non-stressed cultures were used. Samples were enriched directly following inoculation or after storage for seven days at -20 °C. Enrichments were incubated for eight hours and 24 hours. Centrifugation of the sample was evaluated as a pre sample treatment for performance improvement; 160 µl of sample was transferred to the LFA. The result was read after 30 minutes.

## Results

The detection limit of TDH was 125 pg/ml and  $3.3 \times 10^6$  to  $1.9 \times 10^7$  cfu/ml for TDH-positive *V. parahaemolyticus*, strain-dependent. The LFA achieves an inclusivity rate of 81% and exclusivity rate of 100%. For the inclusivity six tdh-gene positive

*V. parahaemolyticus* isolates were tested negative by LFA (**Figure 3**). They are of environmental origin (e.g. seawater and zooplankton) and were tested for TDH production by Latex agglutination test KAP-RPLA (Denka Seiken, Japan). Two of them showed no agglutination (TDH negative) and therefore were excluded in the LFA inclusivity rate calculation.

For fresh food, detection rate of the LFA after 24 hour incubation, in combination with no centrifugation step, was significantly lower than the rate obtained by other methods. In the group of frozen samples, detection after 24 hour enrichment (independent from centrifugation step) was significantly higher than after eight hour enrichment.



**Figure 2:** Test procedure of food experiments. Comparison of LFA method ISO/TS 21872-1:2007

Species	n	Result LFA	Inclusivity/ Exclusivity
<i>V. parahaemolyticus</i>	<i>tdh</i> - gene positive	n = 21 n = 17 positive	Inclusivity 81%
	<i>tdh</i> - gene negative	n = 69 n = 69 negative	Exclusivity 100%
other <i>Vibrio</i> spp.	<i>V. holisae</i> ATCC 33564 <i>V. cholerae</i> Eltor Ogawa NIH 41 <i>V. vulnificus</i> ATCC 27562 <i>V. vulnificus</i> ATCC 33149 <i>V. mimicus</i> ATCC 33653 <i>V. alginolyticus</i> H6533	n = 6 n = 6 negative	Exclusivity 100%
non- <i>Vibrio</i>	<i>E. coli</i> ATCC 25922 <i>C. freundii</i> ATCC 8090 <i>A. hydrophila</i> ATCC 7966 <i>P. shigelloides</i> ATCC 14029	n = 4 n = 4 negative	Exclusivity 100%

Figure 3: Results of inclusivity and exclusivity testing

For both fresh and frozen food types, 100% sensitivity was achieved by LFA after 24 hour enrichment in combination with sample centrifugation. Performance was equivalent to the ISO/TS 21872-1:2007 reference method (100% sensitivity) and time-to-result was achieved four days faster. The preliminary centrifugation treatment of the sample significantly increases the detection rate ( $p=0.035$ ). None of the negative controls were contaminated with TDH-positive *V. parahaemolyticus*, but sporadically with TDH-negative *V. parahaemolyticus*. All negative controls reacted negatively by LFA. Therefore, the specificity of the LFA was 100%.

## Conclusion

Food experiments with artificially contaminated seafood samples showed that TDH-positive *V. parahaemolyticus* was reliably detected in inoculation concentrations of 101 to 102 cfu/g in fresh food and 103 to 104 cfu/g in frozen food after 24 hour incubation (Figure 4).

A GLISA for the detection of pathogenic *V. parahaemolyticus* in food was developed by targeting the toxin TDH. The detection limit of TDH was 125 pg/ml and  $3.3 \times 10^6$  to  $1.9 \times 10^7$  cfu/ml for TDH-positive *V. parahaemolyticus*. In internal studies ( $n=102$ ) a sensitivity of 81% and specificity of 100% was determined for the developed test.

Method of Detection	Type of Food		
	Fresh	Frozen	Total
LFA 8 h/ without centrifugation	<b>100%</b> (94.5% - 100%) <sup>1</sup> n = 53	<b>91.7%</b> (77.5% - 98.2%) <sup>1</sup> n = 36	<b>96.6%</b> (90.4% - 99.3%) <sup>1</sup> n = 89
LFA 8 h/ with centrifugation	<b>100%</b> (94.4% - 100%) <sup>1</sup> n = 52	<b>91.6%</b> (77.5% - 98.2%) <sup>1</sup> n = 36	<b>96.5%</b> (90.4% - 99.3%) <sup>1</sup> n = 88
LFA 24 h/ without centrifugation	<b>85.2%</b> (72.9% - 93.4%) <sup>1</sup> n = 54	<b>100%</b> (92.0% - 100%) <sup>1</sup> n = 36	<b>91.1%</b> (83.2% - 96.1%) <sup>1</sup> n = 90
LFA 24 h/ with centrifugation	<b>100%</b> (94.6% - 100%) <sup>1</sup> n = 54	<b>100%</b> (92.0% - 100%) <sup>1</sup> n = 36	<b>100%</b> (96.7% - 100%) <sup>1</sup> n = 90
Reference method	<b>100%</b> (94.6% - 100%) <sup>1</sup> n = 54	<b>100%</b> (92.0% - 100%) <sup>1</sup> n = 36	<b>100%</b> (96.7% - 100%) <sup>1</sup> n = 90

n number of samples, which were included in statistical analysis  
(viz. excluding samples with invalid test results)

<sup>1</sup> confidence interval

Figure 4: Results of food experiments

Experiments show that *V. parahaemolyticus* in seafood can be detected much faster using lateral flow technology than with traditional methods. Detection was completed in 24 hours with enrichment plus one hour sample pre-treatment and assay performance compared to three to seven days for standard detection methods.

To protect product integrity and improve consumer safety, rapid methods for detection of contamination in food are increasingly important. Tests based on lateral flow technology are user-friendly, does not require specialized, cost-intensive equipment, unlike alternative methods such as PCR. Pathogen detection can therefore be performed more efficiently and cost-effectively.



<sup>1</sup>World Health Organization (2012). Initiative to estimate the Global Burden of Foodborne Diseases. Retrieved March 7, 2013, from [http://www.who.int/foodsafety/foodborne\\_disease/ferg/en/index.html](http://www.who.int/foodsafety/foodborne_disease/ferg/en/index.html)

<sup>2</sup>Centers for Disease Control and Prevention (2012). CDC Estimates of Foodborne Illness in the United States. Retrieved November 20, 2012, from <http://www.cdc.gov/foodborneburden/2011-foodborne-estimates.html>

<sup>3</sup>European Food Safety Authority. Food-borne Zoonotic Diseases. Retrieved November 20, 2012, from <http://www.efsa.europa.eu/en/topics/topic/foodbornezoonoticdiseases.htm>

<sup>4</sup>FAO/WHO [Food and Agriculture Organization of the United Nations/World Health Organization]. 2011. Risk assessment of *Vibrio parahaemolyticus* in seafood: Interpretative summary and Technical report. Microbiological Risk Assessment Series No. 16. Rome. 193pp

<sup>5</sup>Centers for Disease Control and Prevention (2012). CDC Estimates of Foodborne Illness in the United States. Retrieved July 29, 2013, from <http://www.cdc.gov/nczved/divisions/dfbmd/diseases/vibriop/>

# Fed Up with the Challenges of Food Safety?

## Simplify Your Process

### Singlepath® & Duopath® Lateral Flow Tests

#### Convenient, Rapid Food Pathogen Detection

EMD Millipore's Singlepath and Duopath Lateral Flow Tests are immunoassays for detecting pathogens. They offer all the benefits of traditional testing methods with the addition of simplicity, speed, reliability, and convenience. Covering the major pathogens, the tests act as mini-laboratories in the 'pregnancy test' format and always include a built-in control reaction. Use in combination with our granulated culture media to ensure optimal test performance. Singlepath *E. coli* 0157, Singlepath Salmonella, Singlepath Campylobacter, and Duopath Verotoxins rapid tests are AOAC-R1 approved.

**Reliable:** Same accuracy standards as classical detection methods


**Comprehensive:** Tests cover six most relevant pathogens in food

**Fast:** Definite results within 20 - 30 minutes

**Easy-to-use:** Clear yes/no results after simple sample application

**Safe:** Additional positive control and specially adapted enrichment media for precise, reliable results

**Economic:** Rapid results save operating costs by expediting product release while ensuring product reliability



Lateral Flow Tests are available for:

- *E. coli*
- Bacillus cereus Emetic Toxin Marker
- Campylobacter
- L'mono
- Salmonella
- Cereus Enterotoxins
- Legionella

Description	Cat. No.	Each
Singlepath <i>E. coli</i> 0157	EM1.04141.0001	183.08
Singlepath Bacillus cereus Emetic Toxin Marker	EM1.04154.0001	178.00
Singlepath Campylobacter	EM1.04143.0002	348.19
Singlepath L'mono	EM1.04148.0002	225.97
Singlepath Salmonella	EM1.04140.0001	199.00
Duopath Cereus Enterotoxins	EM1.04146.0001	368.23
Duopath Legionella	EM1.04147.0001	305.29

Each test system includes 25 tests.



# Analysis of Free Chlorine in Poultry Chiller Water using DPD-Free Chlorine Reagent

Guidance for improving performance of free chlorine measurements

*Pat Wiese, Applications Chemist, Hach Company*

Poultry processors are required to maintain a free chlorine level sufficient to control pathogenic bacteria found on bird carcasses. Bird carcasses are chilled at 32 to 34°F in a continuous counter-flow chlorinated water system. The challenge to the processor is to maintain a free chlorine level of 0.5 - 5.0 mg/L free chlorine in a chiller water which has a huge and changing chlorine demand load due to fresh carcasses being constantly added. The analysis of free chlorine in these chiller "red" waters is challenging and somewhat controversial as to whether the analysis is accurate. In order to attack this question there are multiple variables that need to be addressed and some assumptions that must be made.

First, the free chlorine concentration is variable throughout the system. The free chlorine concentration will be much higher at the point of chlorination than it will be at the opposite end of the chiller where fresh carcasses are being added. The chlorine concentration will vary at sampling sites between these two points. This needs to be taken into consideration when making chlorine additions or when comparing analysis results between different samples or sampling sites.

Second, the chiller waters are cold and the reaction rate of chlorine with the organic matter present is slow. The rate of reaction depends on the type and structure of the organic matter, temperature, and chlorine concentration. It is important to recognize that the "dynamic equilibrium" existing in the sample may be disrupted or shifted when the sample is removed from the system. A sample removed will slowly warm with the residual free chlorine becoming more reactive. This leads to further reaction of the free chlorine with any existing chlorine demand and hence will lead to decreased free chlorine values if the free chlorine analysis is delayed. This becomes important when samples are transported to a lab, or if samples are being analyzed online by a process analyzer. This can be illustrated by taking a sample and analyzing

for free chlorine. Allow the remaining portion of sample to remain at room temperature for 15 - 30 minutes. Reanalyze the sample. The free chlorine values will usually be lower or entirely gone. This is often interpreted as the free chlorine test not giving reproducible results when, in fact, the free chlorine concentration in the sample has actually changed since the original analysis.

Third, DPD reagents for free chlorine react slowly with other inorganic and organic chlorine species present in the sample. The rate at which they react or "interfere" in the free chlorine determination depends on the form and concentration of these compounds present. The instructions for determining free chlorine using DPD Free Chlorine Reagent are written to read the sample as soon as possible after reagent addition and always to read the sample before one minute. Free chlorine reacts immediately with DPD Chlorine Reagent. The reaction time can be extended to one minute to allow any bubbles to dissipate or for particles present in the sample to settle before taking the reading. This effect of interfering compounds is illustrated when a sample is read at one minute after reagent addition and then reread again after five or ten minutes. The free chlorine value will be higher after five or ten minutes. This is often interpreted as "My free chlorine value keeps increasing!" This continual drifting higher of the free chlorine value should actually be attributed to interference from the other chlorinated compounds present in the sample and the value obtained at one minute or less is considered most accurate.

Fourth, the physical aspects of the sample such as temperature, background red color, particles, and turbidity need to be addressed. The colorimeter or spectrophotometer must always be zeroed on the sample or the sample dilution being analyzed before addition of the DPD reagent. This will automatically subtract the background due to color and turbidity. One precaution should be noted. If bubbles or large particles are present in the sample,



these should be allowed to dissipate or settle out before zeroing the DR 1900 Portable Spectrophotometer, or other Hach Instrument.

Fifth, sampling technique must be considered. Contamination from chlorine demand from previous samples will cause low results. Fats and other organics adhere to container and sample cell walls, reacting with any free chlorine present in the sample. Sampling containers should not be reused unless pretreated for chlorine demand. Glass containers are preferred, but often are not allowed in the processing area. Pre-treat the sample containers to remove any chlorine demand by soaking the container in a dilute bleach solution (1mL commercial bleach to 1L of deionized water) for at least one hour. Rinse each container or cell thoroughly with deionized or organic-free water before using. After testing is completed, rinse the sample cells and sampling containers with deionized water and fill with the weak bleach solution. This will insure that the cells and bottles are ready to use for the next testing period.

These strategies need to be incorporated into all testing protocols in order to develop confidence in the test results and to ensure that the results are useful in making process adjustments and controlling pathogenic bacteria levels. To see the Hach method for determining free chlorine in poultry chiller water visit [vwr.com/hachDR1900](http://vwr.com/hachDR1900).



## DR 1900 Portable Spectrophotometer

It is time to bring your spectrophotometer up to speed with the rest of your life. The new DR 1900 Spectrophotometer from Hach offers you the reliability you've come to know, with additional features that will make it a staple in your lab for the future.

- Small and streamlined footprint that saves valuable bench space and offers true portability
- Large, easy-to-read backlit display
- Intuitive user interface that doesn't require training
- Accepts the widest range of vial sizes
- Over 220 commonly used pre-programmed methods

Description	Cat. No.	Each
DR 1900 Portable Spectrophotometer	10145-568	3250.69



## HQd Water Quality Meters & Probes



Hach's user-friendly meters allow even new users to easily navigate without errors. And only HQd digital probes store calibration history and automatically recognize testing parameters and method settings – reducing possible errors while saving you time and hassle.

The durable design is built to do the job no matter where you put it to work.

Description	Cat. No.	Each
General Purpose pH Kit	10146-094	1982.51
General Purpose pH Kit, High-Performance	10146-096	2111.91
High-Performance General Purpose and LIS	10146-098	2111.91
High-Performance Kit for "Dirty" Samples	10146-100	2163.67
High-Performance Kit for "Challenging" Samples	10146-102	2174.03

## AT1000 Automatic Titrators



Simplify complex titration with the TitraLab AT1000 from Hach. Straight from the box the AT1000 eliminates operator interpretation and manual processes you have come to expect with existing manual titration. The automatic system gives you more reliable results without complicated analysis, all in application specific kits. Visit [vwr.com](http://vwr.com) to see all of the kits available!

Description	Cat. No.	Each
TitraLab AT1000 Titrator, 1 Syringe	10662-252	2390.00
TitraLab AT1000 Titrator, 1 Syringe, 1 Pump	10662-254	2710.00
TitraLab AT1000 Titrator, 1 Syringe, 2 Pumps	10662-256	3030.00
TitraLab AT1000 Titrator, 2 Syringe, 2 Pumps	10662-258	3420.00
TitraLab KF1000, KF, 1 Syringe, 2 Pumps	10662-260	3810.00

## Sentino® Family of Microbiology Systems, Choosing the Right Product is as Easy as 1, 2, 3

The Sentino Family of Microbiology products offers a mix and match selection of products to best suit the economic, ergonomic, and workflow needs in a busy microbiology laboratory. The collections of complimentary products are targeted for evaluating microbial contamination in aqueous samples using MF Technique. Select the items that best fit the needs in your laboratory. Choose disposable filter funnels and our Sentino Pump, or pair the pump with our Sentino Filter Dispenser with individual membrane filters aseptically dispensed at the press of a button. The compact design of the Sentino Microbiology System frees valuable bench-top space and provides flexibility in arranging workspace for optimal efficiency and workflow.

1. Select your device
2. Select your membrane
3. Select your equipment

Description	Cat. No.	Price
Sentino Microbiology Pump	89183-082	Ea./ 1340.61
Sentino Membrane Dispenser	10147-394	Ea./ 1900.00
Sentino Filter Funnel, 100 mL	10147-760	Cs. 100/ 90.08
Sentino Filter Funnel, 250 mL	10147-762	Cs. 80/ 89.67
Sentino Filter Funnel Adaptor	10147-386	Ea./ 141.77
Sentino Filter Funnel, 3-Place Manifold	10147-390	Ea./ 1065.51
Sentino Filter Funnel, 6-Place Manifold	10147-392	Ea./ 1919.32

Visit [vwr.com](http://vwr.com) for a complete listing of all the Sentino products and membranes available.



## Proven Value to Help You Re-Focus Valuable Research Time

With support from **VWRCATALYST**, your business can maintain a competitive edge. Our services re-focus scientific time to help you accelerate productivity and drive scientific innovation.

### Comprehensive services from **VWRCATALYST** include:

- Procurement and Supply Management
- Laboratory and Production Support
- Scientific Support
- Equipment and Instrument Services
- Lean Six Sigma Process Consulting


Over 1,200 **VWRCATALYST** associates are working worldwide today at industry-leading pharmaceutical, biotech, education, industrial, and high-tech production institutions.

### We Enable Science by:

- Powering productivity
- Improving quality, safety, and regulatory compliance
- Reducing total operating costs

Watch an introductory video, download case studies, and more at [vwr.com/vwrcatalyst](http://vwr.com/vwrcatalyst).


# BACTRON Anaerobic Chambers



The Most  
Comfortable  
Chamber

New  
Features

300-900  
Plate Capacity



Model No.	Workspace, cu. ft.	Plate Capacity	Cat. No.	Each
BACTRONEZ	13.7	300	10210-216	11,171.00
BACTRON300	17.6	300	89409-530	24,268.03
BACTRON600	17.6	600	89409-534	28,040.49
BACTRON900	17.6	900	89409-536	29,416.10

Visit [vwr.com](http://vwr.com) for more information.

## SHEL LAB®

## Why Risk Your Company's Reputation?



### Sterileware® Metal Detectable Sampling Tools

- Patented plastic material allows scoop and scoop fragments as small as 1.6mm<sup>3</sup> to be detected with standard metal detection and X-ray equipment
- FDA grade – safe for food, drug and cosmetic contact
- Sterilized by gamma irradiation; individually packaged; lot stamped

Capacity	Cat. No.	Pack of 100
<b>Blue Scoops</b>		
2 oz.	89233-982	182.79
4 oz.	89233-984	273.21
8 oz.	89233-986	462.45
<b>Black Scoops</b>		
2 oz.	89093-476	179.41
4 oz.	89093-480	274.55
8 oz.	89093-478	464.75

**FREE Try-it Offer -  
see page 12 for  
details!**



# TECTA Automated Water Pathogen Detection

**ENDETEC**  **VEOLIA**

**U.S. EPA Approved!**

Even with advances in farming, distribution, processing, and storage practices, *E. coli* outbreaks continue to occur throughout the food chain.

## Stop *E. coli* in Its Tracks

The TECTA Automated Microbiology Platform requires much less manual labor and use of disposable products – the test cartridge is pre-filled with all required test reagents eliminating the need for any handling, dilution, or mixing of reagents for test samples. When results are available in 2-18 hours, appropriate action can be taken if required.

The core TECTA automated microbiology technology can be packaged in a variety of instrument configurations to meet the needs of a wide range of applications – from the

smallest remote municipality to the largest industrial processor.

**Desktop Systems for Municipal and Industrial Quality Assurance** - Samples from municipal water distribution systems, or from points of use in industrial process water systems, are processed by a compact desktop instrument located in close proximity to the point of sampling by personnel who do not require extensive microbiological training.

**Laboratory Automation Systems** - Automate the process of testing large volumes of samples that have been collected in a variety of locations and returned to a central testing laboratory for analysis; eliminating all the sample preparation steps that traditional methods require.



**For more information  
see pages 8-9 in  
this issue!**